



Technology Transition Workshop | *Thomas Hall, Ph.D.*

Overview of the Ibis™ Y-STR, mtDNA and SNP Assays

Assays

- **Ibis™ Y-STR Assay format**
- **Ibis™ mtDNA Assay format**
- **Ibis™ SNP Assay format**

Y-STR Assay Format

Y-STR Markers

- Core minimum haplotype markers + recommended loci DYS437, DYS438 and DYS439




Image adapted from:

<http://www.cstl.nist.gov/biotech/strbase/ystrpos1.htm>


Technology
Transition Workshop 

Additional Y-STR Loci Being Considered

- **DYS456**
- **DYS458**
- **DYS448**
- **DYS635**
- **Y-GATA-H4**

**Included in various tested
assay configurations**


Approach to Y-STR Assay Development



Original Y-STR Assay Layout


Reaction	Primer pair	Locus	Conc (nM)
Multiplex 1	4586	DYS389I	160
	4597	DYS392	160
	4594	DYS391	160
	4602	DYS393	160
	4591	DYS390	160
Multiplex 2	4587	DYS389II-1	200
	4611	DYS438	200
	4615	DYS439	200
4608	DYS437	200	
Single-plex 1	4579	DYS19	250
Single-plex 2	4582	DYS385a/b	250

- 24 samples per plate
- 5-plex, 4-plex and 2 single-plex reactions



Reducing the Size of DYS385a/b

- DYS385a/b has a large product size range
- Shortest primer pair in STRbase has range of 241 – 324 bp
- Our system does best at < 200 bp




- This priming strategy has been tested with 126 samples and confirmed by standard typing

Reducing the Size of *DYS385a/b*

- Two blood-derived DNAs tested at 1 ng each


072109C

DYS385a/b
14: [A93 G46 C5 T13] (157 bp)




SC35495

DYS385a/b
14: [A93 G46 C5 T13] (157 bp)




*Splitting **DYS389I/II** into Separate Products*

- **DYS389** has a duplicated forward primer binding region
- The duplicated primer target spans one repeat region
- The reverse primer and second forward primer span **DYS389I**
- Two products are generated with one primer pair



Splitting *DYS389I/II* into Separate Products


- *DYS389* has a duplicated forward primer binding region
- There is a 4-base difference that can be exploited
- Modified upstream primer won't prime downstream site
- Modified downstream primer won't prime upstream site
- A reverse primer can be made to pair with the upstream forward primer



- This priming strategy has been tested with 126 samples and confirmed by standard typing
- *DYS389II* = *DYS389I* + *DYS389II-1*
 - When AB has *DYS389I* = 13 and *DYS389II* = 29, Ibis™ gets *DYS389II-1* = 16

Splitting *DYS389I/II* into Separate Products

$$DYS389II = 13 + 16 = 29$$



Y-STR Assay Testing


- 95 male population samples from NIST (John Butler) tested
 - 13 loci: 92 samples had previous Y-STR truth data
- 36 male samples at Ibis™ – 16 loci – 34 also typed with Y-Filer
- All loci were concordant

Locus	Number of alleles observed with SNP(s)	Samples tested	%
DYS635	14	36	38.9
DYS389II-1	42	131	32.1
DYS437	28	131	21.4
DYS390	10	131	7.6
DYS385a/b	5	131	3.8
DYS458	1	36	2.8
DYS393	2	131	1.5
DYS438	2	131	1.5
DYS391	1	131	0.8
DYS19	0	131	0.0
DYS389I	0	131	0.0
DYS392	0	131	0.0
DYS439	0	131	0.0
DYS448	0	36	0.0
DYS456	0	36	0.0
Y-GATA-H4	0	36	0.0

Updated Y-STR Assay Layout

Reaction	Primer pair	Locus
Multiplex 1	4586	DYS389I
	4597	DYS392
	4594	DYS391
	4601	DYS393
	4591	DYS390
	4924	DYS458
Multiplex 2	4587	DYS389II-1
	4611	DYS438
	4615	DYS439
	4608	DYS437
	4929	DYS456
Single-plex 1	4579	DYS19
Single-plex 2	4692	DYS385a/b


- Still 24 samples per plate
- **DYS456 and DYS458 have been added**
- **6-plex, 5-plex and 2 single-plex reactions**




Technology
Transition Workshop




6-plex Reaction 1: Sample N31773




5-plex Reaction 2: Sample N31773



Single Reactions 3 and 4: Sample N31773




Single Reactions 3 and 4: Sample N31773



Sensitivity

Template 072109C

- Dilutions from 1 ng down to 7.8 pg per reaction produced full profiles down to 62.5 pg/reaction (250 pg for four reactions/sample)





Template Quantity (pg)	DYS19	DYS385a/b	DYS389I	DYS389II-1	DYS390	DYS391	DYS392	DYS393	DYS437	DYS438	DYS439	DYS456	DYS458
1000	14	13, 14	12	17	23	10	11	14	16	10	11	14	16
	14	13, 14	12	17	23	10	11	14	16	10	11	14	16
500	14	13, 14	12	17	23	10	11	14	16	10	11	14	16
	14	13, 14	12	17	23	10	11	14	16	10	11	14	16
250	14	13, 14	12	17	23	10	11	14	16	10	11	14	16
	14	13, 14	12	17	23	10	11	14	16	10	11	14	16
125	14	13, 14	12	17	23	10	11	14	16	10	11	14	16
	14	13, 14	12	17	23	10	11	14	16	10	11	14	16
	14	13, 14	12	17	23	10	11	14	16	10	11	14	16
	14	13, 14	12	17	23	10	11	14	16	10	11	14	16
62.5 pg	14	13, 14	12	17	23	---	11	12, 14	16	10	11	14	16
	14	13, 14	12	17	23	10	11	12, 14	14, 16	10	11	14	16
	14	13, 14	12	17	23	10	11	14	16	10	11	14	16
15.6	14	13, 14	12	17	23	10	11	14	16	10	11	14	16
	14	13, 14	12	17	23	10	11	14	16	10	11	13, 14, 15	16
7.8	14	13, 14	---	16, 17	23	10	11	14	16	10	11	15	16
	14	13, 14	---	16, 17	23	10	---	14	---	10	---	14	16
Negative	---	---	---	---	---	---	---	---	---	---	---	---	---
	---	---	---	---	---	---	---	---	---	---	---	---	---
	---	---	---	---	---	---	---	---	---	---	---	---	---
	---	---	---	---	---	---	---	---	---	---	---	---	---
	---	---	---	---	---	---	---	---	---	---	---	---	---
	---	---	---	---	---	---	---	---	---	---	---	---	---

Species Specificity

- **Human blood-derived DNA sample was tested in triplicate in the presence of a 10-fold excess of exogenous DNA**
- **1 ng of human DNA per reaction**
- **10 ng exogenous DNA**
 - **Dog (male American Eskimo – buccal swab)**
 - **Cat (male long-hair, buccal swab)**
 - **Candida albicans (yeast)**
 - **Aspergillus oryzae (environmental filamentous fungus)**
 - **Escherichia coli (gram negative bacterium)**
 - **Staphylococcus aureus (gram positive bacterium)**
- **All tests with exogenous DNA gave a full profile**


Species Specificity – 6-plex Reaction 1

1 ng human DNA + 10 ng dog DNA




Species Specificity – 5-plex Reaction 2

1 ng human DNA + 10 ng dog DNA







10 ng dog DNA



Species Specificity – Single Reactions 3 and 4

1 ng human DNA + 10 ng dog DNA



Species Specificity

Sample	DYS19	DYS385a/b	DYS389I	DYS389II-1	DYS390	DYS391	DYS392	DYS393	DYS437	DYS438	DYS439	DYS456	DYS458
Dog + Human DNA	16	14, 15	13	18	24	11	11	13	15	10	12	15	18
	16	14, 15	13	18	24	11	11	13	15	10	12	15	18
	16	14, 15	13	18	24	11	11	13	15	10	12	15	18
Cat + Human DNA	16	14, 15	13	18	24	11	11	13	15	10	12	15	18
	16	14, 15	13	18	24	11	11	13	15	10	12	15	18
	16	14, 15	13	18	24	11	11	13	15	10	12	15	18
Staphylococcus aureus + Human DNA	16	14, 15	13	18	24	11	11	13	15	10	12	15	18
	16	14, 15	13	18	24	11	11	13	15	10	12	15	18
	16	14, 15	13	18	24	11	11	13	15	10	12	15	18
Escherichia coli + Human DNA	16	14, 15	13	18	24	11	11	13	15	10	12	15	18
	16	14, 15	13	18	24	11	11	13	15	10	12	15	18
	16	14, 15	13	18	24	11	11	13	15	10	12	15	18
Candida albicans + Human DNA	16	14, 15	13	18	24	11	11	13	15	10	12	15	18
	16	14, 15	13	18	24	11	11	13	15	10	12	15	18
	16	14, 15	13	18	24	11	11	13	15	10	12	15	18
Aspergillus oryzae + Human DNA	16	14, 15	13	18	24	11	11	13	15	10	12	15	18
	16	14, 15	13	18	24	11	11	13	15	10	12	15	18
	16	14, 15	13	18	24	11	11	13	15	10	12	15	18
Human DNA alone	16	14, 15	13	18	24	11	11	13	15	10	12	15	18
	16	14, 15	13	18	24	11	11	13	15	10	12	15	18

Exogenous DNAs did not interfere with full profile detection

Cat alone (10 ng/reaction)	---	---	---	---	---	---	---	---	---	---	---	---	---
	---	---	---	---	---	---	---	---	---	---	---	---	---
	---	---	---	---	---	---	---	---	---	---	---	---	---
Staphylococcus aureus alone (10 ng/reaction)	---	---	---	---	---	---	---	---	---	---	---	---	---
	---	---	---	---	---	---	---	---	---	---	---	---	---
	---	---	---	---	---	---	---	---	---	---	---	---	---
Escherichia coli alone (10 ng/reaction)	---	---	---	---	---	---	---	---	---	---	---	---	---
	---	---	---	---	---	---	---	---	---	---	---	---	---
	---	---	---	---	---	---	---	---	---	---	---	---	---
Candida albicans alone (10 ng/reaction)	---	---	---	---	---	---	---	---	---	---	---	---	---
	---	---	---	---	---	---	---	---	---	---	---	---	---
	---	---	---	---	---	---	---	---	---	---	---	---	---
Aspergillus oryzae alone (10 ng/reaction)	---	---	---	---	---	---	---	---	---	---	---	---	---
	---	---	---	---	---	---	---	---	---	---	---	---	---
	---	---	---	---	---	---	---	---	---	---	---	---	---
Water	---	---	---	---	---	---	---	---	---	---	---	---	---
	---	---	---	---	---	---	---	---	---	---	---	---	---
	---	---	---	---	---	---	---	---	---	---	---	---	---

Software and Database in Place for Y-STRs

The screenshot displays the IbisTrack software interface for Y-STR analysis. The main window shows a plate analysis for Plate P05014932, Sample 2, dated 05-Apr-2009. The plate grid at the top lists various Y-STR markers across 12 wells. The main analysis area features a bar chart for each marker, with the following observed allele sizes (D=I):

Marker	Observed Allele Size (D=I)
DYS438	12
DYS439	12
DYS385a/b	12, 14
DYS19	14
DYS389I	14
DYS389II-1	17
DYS390	24
DYS391	9
DYS392	13
DYS393	13, 14
DYS437	15

The detailed view for Well 38 (D02) shows the following data:

- Well: 38 (D02), PP 4692 (DYS385-A-B_AC022496_29490_2)
- Locus: **DYS385a/b**
- Allele: **12**
- Base count: A87 G44 C5 T13
- Top strand: Abundance: 4368.2752, Predicted mass: 47407.9023, Observed mass: 47408.8683, Match error: -20.377 ppm
- Bottom strand: Abundance: 4802.5964, Predicted mass: 44844.7343, Observed mass: 44846.1918, Match error: -32.502 ppm

A table below the detailed view shows the genotype for other markers in the well:

Locus	Allele	Allele
DYS19	14	
DYS385a/b	12	14
DYS389I	14	
DYS389II-1	17	
DYS390	24	
DYS391	9	
DYS392	13	
DYS393	13	14 (T→C)
DYS437	15	
DYS438	12	
DYS439	12	

The interface also includes a 'Comments for analysis report' section at the bottom.


mtDNA Assay Format

Outline


- **Ibis™ mtDNA Assay format**
- **Base composition analysis of mtDNA samples**
 - Data processing and analysis
- **Information content relative to sequencing**
- **Heteroplasmy detection**
- **Sensitivity**
- **Reproducibility**

mtDNA Tiling Assay Format

- 24 primer pairs cover amplified coordinates HV1 15924..16428 and HV2 31..576
- Target most highly-conserved positions on 3' ends of primers
- Grouped into eight triplexed sets (colored grouping) by maximum spatial separation and suitable mass separation of products



mtDNA Tiling Assay Format




- Primers grouped to maximize target site separation
- PCR reactions performed with short extension cycle (five seconds)
- Product masses resolve from each other in triplex groupings
- Product sizes range from 85 to 140 bp; all but three are <150 bp
- Relative primer pair concentrations in triplexes have been adjusted to favor simultaneous amplification of all products

Base Composition Analysis


- 1. PCR**
- 2. Desalting**
- 3. ESI-TOF mass spectrometry**
- 4. Raw spectrum processing/deconvolution**
- 5. Base composition assignment/profile development**

PCR: Three Primer Pairs Per Reaction




Desalting of PCR Reactions

Magnetic bead anion exchange



ESI-TOF Mass Spectrometry

- Three primer pairs per reaction
- Complex raw spectrum




Technology
Transition Workshop



ESI-TOF Mass Spectrometry

- Three primer pairs per reaction
- Complex raw spectrum



- Multiple peaks for each mass
- 3 (or more) ds products per spectrum
- 6 DNA strands per spectrum


Noise baseline

Technology
Transition Workshop



ESI-TOF Mass Spectrometry

- Three primer pairs per reaction
- Complex raw spectrum



- Multiple peaks for each mass
- 3 (or more) ds products per spectrum
- 6 DNA strands per spectrum


Baseline corrected spectrum used for deconvolution

Technology
Transition Workshop



Raw Spectrum Processing


- Internal mass standards bracket the spectrum for accurate calibration of the measurements before deconvolution



Technology
Transition Workshop



Spectral Deconvolution




Technology
Transition Workshop




Deconvolution to Masses

- Deconvolution results in one final measurement per molecular species




DS DNA Strand Association

- Forward and reverse strands of a double stranded DNA can be associated by mass




DS DNA Strand Association

- Forward and reverse strands of a double stranded DNA can be associated by mass




Base Composition Assignment



[A40 G9 C40 T19] [A47 G18 C25 T30] [A49 G17 C31 T40]





Base Composition Assignment



Final Product Assignment





Final Product Assignment





15893..16012
[A47 G18 C25 T30]

Final Product Assignment



15893..16012
[A47 G18 C25 T30]



Final Product Assignment



15893..16012
[A47 G18 C25 T30]

16231..16338
[A40 G9 C40 T19]

Final Product Assignment





15893..16012
[A47 G18 C25 T30]

16231..16338
[A40 G9 C40 T19]

Technology
Transition Workshop



Final Product Assignment



15893..16012
[A47 G18 C25 T30]


16231..16338
[A40 G9 C40 T19]

154..290
[A49 G17 C31 T40]


Technology
Transition Workshop




Eight Spectra Per Sample





Coverage Map Development



P20000040-AF-12NQ composite



Base Composition Profile




Base composition profile


Base composition
rCRS coordinates
Primer pair

2901:	15893..16012:	A47	G18	C25	T30
2925:	15937..16041:	A35	G14	C24	T32
2899:	15985..16073:	A26	G15	C21	T27
2898:	16025..16119:	A26	G17	C26	T26
2897:	16055..16155:	A31	G13	C30	T27
2896:	16102..16224:	A45	G13	C42	T23
2895:	16130..16224:	A36	G7	C33	T19
2893:	16154..16268:	A44	G7	C46	T18
2892:	16231..16338:	A40	G9	C40	T19
2891:	16256..16366:	A37	G9	C41	T24
2890:	16318..16402:	A20	G14	C30	T21
2889:	16357..16451:	A21	G17	C36	T21
2902:	5..97:	A19	G24	C24	T26
2903:	20..139:	A24	G34	C29	T33
2904:	83..187:	A23	G21	C29	T32
2905:	113..245:	A39	G18	C28	T48
2906:	154..290:	A49	G17	C31	T40
2908:	204..330:	A42	G16	C35	T32
2908:	204..330:	A42	G16	C36	T32
2908:	204..330:	A42	G16	C37	T32
2908:	204..330:	A42	G16	C38	T32
2907:	239..363:	A43	G11	C46	T23
2907:	239..363:	A43	G11	C47	T23
2907:	239..363:	A43	G11	C48	T23
2907:	239..363:	A43	G11	C49	T23
2923:	262..390:	A47	G10	C50	T20
2923:	262..390:	A47	G10	C51	T20
2923:	262..390:	A47	G10	C52	T20
2923:	262..390:	A47	G10	C53	T20
2910:	331..425:	A33	G9	C27	T26
2918:	367..463:	A27	G8	C32	T30
2912:	409..521:	A32	G7	C48	T26
2913:	464..603:	A44	G10	C63	T23


Multiple Assays Per Plate




Multiple Mass Assignments




Multiple Samples Per Plate




Multiple Profiles Per Plate



2901: 15893..16012: A47 G18 C25 T30
 2925: 15937..16041: A35 G14 C24 T32
 2899: 15985..16073: A26 G15 C21 T27
 2898: 16025..16119: A26 G17 C26 T26
 2897: 16055..16155: A31 G13 C30 T27
 2896: 16102..16224: A45 G13 C42 T23
 2895: 16130..16224: A36 G7 C33 T19
 2893: 16154..16268: A44 G7 C46 T18
 2892: 16231..16338: A40 G9 C40 T19
 2891: 16256..16366: A37 G9 C41 T24
 2890: 16318..16402: A20 G14 C30 T21
 2889: 16357..16451: A21 G17 C36 T21
 2902: 5..97: A19 G24 C24 T26
 2903: 20..139: A24 G34 C29 T33
 2904: 83..187: A23 G21 C29 T32
 2905: 113..245: A39 G18 C28 T48
 2906: 154..290: A49 G17 C31 T40
 2908: 204..330: A42 G16 C35 T32
 2908: 204..330: A42 G16 C36 T32
 2908: 204..330: A42 G16 C37 T32
 2908: 204..330: A42 G16 C38 T32
 2907: 239..363: A43 G11 C46 T23
 2907: 239..363: A43 G11 C47 T23
 2907: 239..363: A43 G11 C48 T23
 2907: 239..363: A43 G11 C49 T23
 2923: 262..390: A47 G10 C50 T20
 2923: 262..390: A47 G10 C51 T20
 2923: 262..390: A47 G10 C52 T20
 2923: 262..390: A47 G10 C53 T20
 2910: 331..425: A33 G9 C27 T26
 2916: 367..463: A27 G8 C32 T30
 2912: 409..521: A32 G7 C48 T26
 2913: 464..603: A44 G10 C63 T23



2901: 15893..16012: A46 G19 C25 T30
 2925: 15937..16041: A35 G14 C24 T32
 2899: 15985..16073: A26 G15 C21 T27
 2898: 16025..16119: A26 G17 C28 T24
 2897: 16055..16155: A32 G12 C30 T27
 2896: 16102..16224: A46 G12 C41 T24
 2895: 16130..16224: A36 G7 C33 T19
 2893: 16154..16268: A44 G7 C45 T19
 2892: 16231..16338: A40 G9 C40 T19
 2891: 16256..16366: A37 G9 C41 T24
 2890: 16318..16402: A20 G14 C30 T21
 2889: 16357..16451: A22 G16 C36 T21
 2902: 5..97: A19 G24 C24 T26
 2903: 20..139: A24 G34 C29 T33
 2904: 83..187: A23 G21 C30 T31
 2905: 113..245: A39 G18 C31 T45
 2906: 154..290: A48 G18 C35 T36
 2908: 204..330: A42 G16 C40 T30
 2907: 239..363: A43 G11 C49 T23
 2923: 262..390: A47 G10 C53 T20
 2910: 331..425: A33 G9 C27 T26
 2916: 367..463: A27 G8 C32 T30
 2912: 409..521: A32 G7 C48 T26
 2913: 464..603: A45 G10 C68 T23
 2913: 464..603: A45 G10 C69 T23



2901: 15893..16012: A47 G18 C25 T30
 2925: 15937..16041: A35 G14 C24 T32
 2899: 15985..16073: A26 G15 C21 T27
 2898: 16025..16119: A26 G17 C27 T25
 2897: 16055..16155: A32 G12 C29 T28
 2897: 16055..16155: A31 G13 C29 T28
 2896: 16102..16224: A46 G12 C42 T23
 2896: 16102..16224: A45 G13 C42 T23
 2895: 16130..16224: A36 G7 C33 T19
 2893: 16154..16268: A44 G7 C46 T18
 2892: 16231..16338: A39 G10 C40 T19
 2891: 16256..16366: A36 G10 C42 T23
 2890: 16318..16402: A20 G14 C30 T21
 2889: 16357..16451: A21 G17 C36 T21
 2902: 5..97: A20 G23 C24 T26
 2903: 20..139: A25 G33 C29 T33
 2904: 83..187: A23 G21 C29 T32
 2905: 113..245: A39 G18 C29 T47
 2906: 154..290: A48 G18 C32 T39
 2908: 204..330: A42 G16 C39 T32
 2907: 239..363: A43 G11 C50 T23
 2907: 239..363: A43 G11 C51 T23
 2923: 262..390: A47 G10 C54 T20
 2923: 262..390: A47 G10 C55 T20
 2910: 331..425: A33 G9 C27 T26
 2916: 367..463: A27 G8 C32 T30
 2912: 409..521: A32 G7 C48 T26
 2913: 464..603: A44 G10 C63 T23


Technology
 Transition Workshop



Tiling Compared to Sequencing

- **1266 unique tiling region sequences were selected from GenBank genomes**
 - Each sequence differed from all others by at least one base
 - C-stretch length differences ignored
- **Sequences converted to tiling base compositions**
- **Cross-compared for minimum differences using mtDNA search algorithm**
 - Ignores C-stretch length differences
 - Corrects for primer pair overlaps
- **94.2% of unique tiling region sequences were uniquely discriminated by tiling assay**

Tiling Compared to Sequencing




- Tiling assay amplifies 15924-16428 and 31-576
- Minimum HVI +HVII sequences 16024-16365 and 73-340

For the same set of 1266 unique sequences spanning mtDNA tiling coordinates:

- 94.2% can be differentiated with the tiling assay
- 90.2% can be differentiated by sequencing HVI 16024-16365 and HV2 73-340


Length Heteroplasmy Detection

Example = USA.FBI.000009, region 16102..16224 *



Good sequence up to C stretch, then multiple overlapping reads


Opposite read



MS detects multiple products while sequencing requires a pure sample for maximum information content

* Happened to be same sample given as a blinded sample by Bruce Budowle in 2003. This was suggested by our assay in the FBI validation, then confirmed.

Length Heteroplasmy Example



- 2901: 15893..16012: A47 G18 C25 T30
- 2925: 15937..16041: A35 G14 C24 T32
- 2899: 15985..16073: A25 G16 C21 T27
- 2898: 16025..16119: A25 G18 C27 T25
- 2897: 16055..16155: A31 G12 C30 T28
- 2896: 16102..16224: A44 G12 C42 T24
- 2896: 16102..16224: A44 G12 C43 T24
- 2896: 16102..16224: A44 G12 C44 T24
- 2896: 16102..16224: A44 G12 C45 T24
- 2895: 16130..16224: A35 G7 C33 T19
- 2895: 16130..16224: A35 G7 C34 T19
- 2895: 16130..16224: A35 G7 C35 T19
- 2895: 16130..16224: A35 G7 C36 T19
- 2893: 16154..16268: A43 G7 C47 T17
- 2893: 16154..16268: A43 G7 C48 T17
- 2893: 16154..16268: A43 G7 C49 T17
- 2893: 16154..16268: A43 G7 C50 T17
- 2892: 16231..16338: A40 G9 C39 T20
- 2891: 16256..16366: A37 G9 C41 T24
- 2890: 16318..16402: A20 G14 C31 T20
- 2889: 16357..16451: A21 G17 C36 T21
- 2902: 5..97: A19 G24 C24 T26
- 2903: 20..139: A24 G34 C29 T33
- 2904: 83..187: A23 G21 C30 T31
- 2905: 113..245: A39 G18 C30 T46
- 2905: 113..245: A39 G18 C29 T47
- 2906: 154..290: A48 G18 C32 T39
- 2906: 154..290: A48 G18 C31 T40
- 2908: 204..330: A42 G16 C39 T32
- 2908: 204..330: A42 G16 C40 T32
- 2907: 239..363: A43 G11 C49 T24
- 2907: 239..363: A43 G11 C50 T24
- 2923: 262..390: A47 G10 C53 T21
- 2923: 262..390: A47 G10 C54 T21
- 2910: 331..425: A33 G9 C27 T26
- 2916: 367..463: A27 G8 C32 T30
- 2912: 409..521: A32 G7 C48 T26
- 2913: 464..603: A43 G11 C63 T23

Technology
Transition Workshop



SNP Heteroplasmy Detection

From sequence profile

AF-4: 16024-16365

C 16176 N

T 16362 C

C/T heteroplasmy


C/T heteroplasmy

Observed profile

16102..16224: A45 G13 C41 T24
 16102..16224: A45 G13 C40 T25
 16130..16224: A36 G7 C33 T19
 16130..16224: A36 G7 C32 T20

Calculated from truth key


16102..16224: A45 G13 C40 T24 N
 16130..16224: A36 G7 C32 T19 N



From blinded sample AF-4

Sensitivity

- Analysis of five templates in dilution-to-extinction
- Sensitivity ranged from 0.1 pg to 6.3 pg template/primer pair*
- All templates had full profile at 6.3 pg or below per reaction
- Sensitivity criteria for standard QC plate set at 25 pg/reaction



* Based upon total nuclear DNA quantification with Quantifiler®


Technology Transition Workshop 

Reproducibility

- 3,331 trials of one positive control template run over the course of 27 months
- 500 pg template/reaction
- Automated data analysis
- 79,944 expected ds assignments
- 3,298 (99%) full, correct profiles

Distribution of mass measurement deviations for 159,688 DNA strand assignments

Ave. error magnitude was 10.12 ± 8.04 ppm



Summary


- **1051 nucleotide positions covered by 24 primer pairs**
- **Accurate mass measurements and biochemical strategy allow mitochondrial base composition profiles to be developed**
- **Discrimination power is about 94% that of sequencing same region**
- **Discrimination power over tiling region can be greater than sequencing over minimum HV1 and HV2 ranges**
- **Base composition profiles can be compared to each other and to sequence profiles**
- **Databases can be searched and subjected to same type of statistics as a sequence database**
- **Mass spectrometry can resolve heteroplasmy (mixtures)**

SNP Assay Format


Objective

- **PCR/ESI-MS-based assay for human autosomal SNP analysis**
- **Exclude non-contributors to a DNA sample**
- **A random profile match should have very low probability**
- **Minimize population bias**
 - **Multiple markers with about 50% heterozygosity**
 - **Low F_{st} (distribution same in all populations)**
 - **Low detectable genetic linkage (low linkage disequilibrium)**
 - **Use product rule for probability estimates**
 - **Use global q or no q correction for population substructure**
- **40 independent markers**


Ideal Bi-allelic Markers

One from father  One from mother 

Each autosomal chromosome comes as a pair




Bi-allelic markers have one of two states



Any marker for any individual can be both 'A' or both 'B' or one of each

Ideal Bi-allelic Markers

There's a 50% chance of getting either allele on the father's chromosome



And a 50% chance of getting either allele on the mother's chromosome


Each marker is like a coin toss

Coin images ©Microsoft

Technology
Transition Workshop 

Ideal Bi-allelic Markers

There's a 50% chance of getting either allele on the father's chromosome




And a 50% chance of getting either allele on the mother's chromosome

Each marker is like a coin toss


That means $50\% \times 50\% = 25\%$ chance of getting two 'heads' (A allele)

Coin images ©Microsoft

Technology
Transition Workshop 

Ideal Bi-allelic Markers

There's a 50% chance
of getting either
allele on the father's
chromosome



And a 50% chance of
getting either allele
on the mother's
chromosome


Each marker is like a
coin toss

And $50\% \times 50\% = 25\%$ chance of getting two 'tails' (B allele)

Coin images ©Microsoft

Technology
Transition Workshop 

Ideal Bi-allelic Markers




Each marker is like a
coin toss
And $50\% \times 50\% = 25\%$ chance of getting A + B

Coin images ©Microsoft

Technology
Transition Workshop 

Ideal Bi-allelic Markers



Each marker is like a
coin toss

And $50\% \times 50\% = 25\%$ chance of getting A + B

Plus $50\% \times 50\% = 25\%$ chance of getting B + A

= 50% chance of being heterozygous

Coin images ©Microsoft


Technology
Transition Workshop




Ideal Bi-allelic Markers

For any two individuals,
the random odds they will match at any one locus


is




Ideal Bi-allelic Markers



Random match probability for one marker
0.375



Random match probability for another marker
0.375




Random match probability for two-marker profile

0.141

Technology Transition Workshop



Ideal Bi-allelic Markers



Random match probability for four-marker profile

$$= 0.375 \times 0.375 \times 0.375 \times 0.375$$

$$= 0.375^4$$

$$= 0.0198$$

Close to 99% with four perfectly distributed bi-allelic markers

Ideal Bi-allelic Markers

- **With 40 unlinked and perfectly-distributed SNPs, the random match probability would ideally be 0.37540, or 9.15×10^{-18}**

Kidd-40 SNP Panel

- Large collection of SNP positions with data for three major population groups
- Subpanel identified with low F_{st} and high heterozygosity
- Evaluated subpanel over seven populations
- 73 SNPs with $F_{st} < 0.02$ over seven populations
- 40 final SNPs with $F_{st} < 0.06$ and heterozygosity > 0.4 across 40 populations around the world
 - Reference: Pakstis, A., Speed, W., Kidd, J., Kidd, K. “Candidate SNPs for a Universal Individual Identification Panel.” *Human Genetics* 121 (3 – 4) (May 2007): 305 – 317.


The Kidd Approach

90,483 SNPs
from ABI
catalog


Frequency data
for European,
African and
Chinese/Japanese

The Kidd Approach

90,483 SNPs
from ABI
catalog




Informatically
selected




436 SNPs




Tested in
TaqMan[®] SNP
assays across 7
populations



73 SNPs




40 SNPs



Kidd Population Coverage

2070 total individuals
40 total populations
Several major global regions




Global $F_{st} \leq 0.06$ all populations
Average heterozygosity > 0.4
Median LD = 0.010 (ave. = 0.029)

Technology
Transition Workshop



Resolving Power of Kidd SNP Panel

- 'Perfect' 40-SNP panel: random match with probability of 9.15×10^{-18}
- Kidd panel has ave. of about 1×10^{-15} match probability
- Within populations, ave. match probability ranged from 10^{-12} to 10^{-16}



- 1,568 full profiles
- 40 populations
- All pairwise profile-to-profile comparisons
- 1,228,528 pairwise comparisons
- Most people will differ at 15 – 16 loci

Why Use Mass Spectrometry?

- **Unified platform for major DNA forensics applications**
 - Mitochondrial DNA profiling
 - STR analysis
 - Autosomal SNP analysis
 - SNPs, STRs and/or mtDNA can be analyzed automatically on one instrument in the same run
- **High degree of accuracy**
 - Potential rare variant will be resolved rather than missed (could a C \leftrightarrow T position ever present an 'A' or a 'G'?)

Initial Development: Eight 5-plex Reactions

	1	2	3	4	5	6	7	8	9	10	11	12
A		UT57318	WT51362	WA29594	JT51471	OT05897	PT84223	PT84232	GT37778	GT37900	TT51422	ZT80786
B	UT57300	WT51342	WT51373	WA29612	JT51499	OT05898	PT84224	PT84234	GT37812	GT37913	TT51435	ZT80815
C	UT57301	WT51343	WT51378	ZT81387	OT05888	OT05899	PT84225	PT84236	GT37828	JT52076	TT51483	ZT80826
D	UT57302	WT51345	WT51381	MT94859	OT05890	OT05901	PT84226	PT84239	ZT80932	OT07280	TT51511	ZT80863
E	UT57303	WT51354	WT51386	MT94866	OT05892	PT84214	PT84227	PT84240	GT37862	PT85612	TT51530	ZT80865
F	UT57310	WT51355	BC11352	MT94868	OT05893	PT84215	PT84228	PT84241	GT37864	PT85658	ZT80731	ZT80869
G	UT57312	WT51358	MT97172	MT94869	OT05894	PT84216	PT84230	PT84242	GT37869	TT51399	ZT80737	ZT80870
H	UT57317	WT51359	WA29584	MT94875	OT05896	PT84222	PT84231	PT84243	GT37888	TT51407	ZT80782	ZT80925


Add 5 μ L template to each well of a plate and thermocycle

Sample 1
Sample 2
Sample 3
Sample 4
Sample 5
Sample 6
Sample 7
Sample 8
Sample 9
Sample 10
 Positive
 Negative





	1	2	3	4	5	6	7	8	9	10	11	12
A	RS1058083 + RS13182883 + RS1821380 + RS214955 + RS7704770											
B	RS10488710 + RS1478829 + RS7205345 + RS985492 + RS987640											
C	RS1410059 + RS1554472 + RS2073383 + RS279844 + RS6444724											
D	RS13134862 + RS1523537 + RS560681 + RS6811238 + RS9951171											
E	RS2272998 + RS321198 + RS3780962 + RS445251 + RS6591147											
F	RS1019029 + RS12997453 + RS1358856 + RS2503107 + RS740598											
G	RS1109037 + RS2567608 + RS315791 + RS447818 + RS7229946											
H	RS10092491 + RS13218440 + RS1336071 + RS338882 + RS7520386											

Each sample is distributed across one column of an assay plate





Initial Multiplexing Results







Initial Multiplexing Results




Initial Multiplexing Results




Initial Multiplexing Results




Initial Intra-locus Strand Balance




Adjusted Multiplexing Results




Adjusted Multiplexing Results




Adjusted Multiplexing Results



Adjusted Multiplexing Results



Improved Intra-locus Strand Balance




Sensitivity

Full profile
obtained at
63 pg/reaction

Reaction quality
degraded at
125 pg and
below

Locus	2000	1000	500	250	125	63	31	16	8
RS10092491	C, T	C, T	C, T	C, T	C, T	C, T	C, T	C, T	no data
RS1019029	C, T	C, T	C, T	C, T	C, T	C, T	C, T	C, ---	C, T
RS10488710	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---
RS1058083	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G
RS1109037	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, ---	A, ---
RS12997453	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---
RS13134862	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, ---	no data
RS13182883	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---	A, G	G, ---
RS13218440	A, ---	A, ---	A, ---	A, ---	A, ---	A, ---	A, ---	A, ---	A, ---
RS1336071	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, ---	A, ---
RS1358856	A, C	A, C	A, C	A, C	A, C	A, C	A, C	no data	C, ---
RS1410059	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---
RS1478829	A, T	A, T	A, T	A, T	A, T	A, T	A, ---	T, ---	A, T
RS1523537	C, T	C, T	C, T	C, T	C, T	C, T	C, T	C, T	no data
RS1554472	C, T	C, T	C, T	C, T	C, T	C, T	C, T	C, T	no data
RS1821380	C, G	C, G	C, G	C, G	C, G	C, G	C, G	G, ---	C, G
RS2073383	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---
RS214955	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, ---
RS2272998	C, G	C, G	C, G	C, G	C, G	C, G	C, G	C, G	C, ---
RS2503107	A, C	A, C	A, C	A, C	A, C	A, C	A, C	A, C	C, ---
RS2567608	A, G	A, G	A, G	A, G	A, G	A, G	A, ---	A, G	no data
RS279844	A, ---	A, ---	A, ---	A, ---	A, ---	A, ---	A, ---	A, ---	A, ---
RS315791	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---
RS321198	C, T	C, T	C, T	C, T	C, T	C, T	C, T	C, T	no data
RS338882	T, ---	T, ---	T, ---	T, ---	T, ---	T, ---	T, ---	T, ---	T, ---
RS3780962	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---	no data	C, ---
RS445251	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---
RS447818	A, G	A, G	A, G	A, G	A, G	A, G	G, ---	G, ---	A, G
RS560681	A, ---	A, ---	A, ---	A, ---	A, ---	A, ---	A, ---	A, ---	A, ---
RS6444724	C, T	C, T	C, T	C, T	C, T	C, T	C, ---	C, T	no data
RS6591147	C, T	C, T	C, T	C, T	C, T	C, T	C, ---	C, T	C, T
RS6811238	T, ---	T, ---	T, ---	T, ---	T, ---	T, ---	T, ---	T, ---	T, ---
RS7205345	C, G	C, G	C, G	C, G	C, G	C, G	C, G	C, G	C, G
RS7229946	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---
RS740598	A, ---	A, ---	A, ---	A, ---	A, ---	A, ---	A, ---	A, ---	no data
RS7520386	A, G	A, G	A, G	A, G	A, G	A, G	A, ---	A, G	A, G
RS7704770	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---
RS985492	T, ---	T, ---	T, ---	T, ---	T, ---	T, ---	T, ---	T, ---	T, ---
RS987640	T, ---	T, ---	T, ---	T, ---	T, ---	T, ---	T, ---	T, ---	T, ---
RS9951171	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, ---	A, G





Inter-Locus 5-plex Balance in Dilutions



Species Specificity

- Human blood-derived DNA sample was tested in duplicate in the presence of 10-fold excess of exogenous DNA
- 1 ng of human DNA per reaction
- 10 ng exogenous DNA
 - Dog (male American Eskimo – buccal swab)
 - Cat (male long-hair, buccal swab)
 - Candida albicans (yeast)
 - Aspergillus oryzae (environmental filamentous fungus)
 - Escherichia coli (gram negative bacterium)
 - Staphylococcus aureus (gram positive bacterium)
- All tests with exogenous DNA gave a full profile





Species Specificity



Species Specificity

With dog DNA





With *Aspergillus oryzae* DNA



Species Specificity

With dog DNA





With *Aspergillus oryzae* DNA



Species Specificity

With dog DNA

With *Aspergillus oryzae* DNA



Species Specificity

10-fold excess of bacterial DNA

Pos

Pos

E. coli

E. coli

E. coli + human

E. coli + human

S. aureus


S. aureus

S. aureus + human

S. aureus + human

neg

neg



Species Specificity

Full profiles were obtained for all replicates in the presence of 10-fold excess of exogenous DNA from six different sources



Locus	Dog	Cat	C. alb	A. ory	E. coli	S. aur	Dog	Cat	C. alb	A. ory	E. coli	S. aur
RS10092491	C, T	C, T	C, T	C, T	C, T	C, T	C, T	C, T	C, T	C, T	C, T	C, T
RS1019029	C, T	C, T	C, T	C, T	C, T	C, T	C, T	C, T	C, T	C, T	C, T	C, T
RS10488710	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---
RS1058083	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G
RS1109037	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G
RS12997453	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---
RS13134862	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G
RS13182883	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---
RS13218440	A, ---	A, ---	A, ---	A, ---	A, ---	A, ---	A, ---	A, ---	A, ---	A, ---	A, ---	A, ---
RS1336071	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G
RS1358856	A, C	A, C	A, C	A, C	A, C	A, C	A, C	A, C	A, C	A, C	A, C	A, C
RS1410059	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---
RS1478829	A, T	A, T	A, T	A, T	A, T	A, T	A, T	A, T	A, T	A, T	A, T	A, T
RS1523537	C, T	C, T	C, T	C, T	C, T	C, T	C, T	C, T	C, T	C, T	C, T	C, T
RS1554472	C, T	C, T	C, T	C, T	C, T	C, T	C, T	C, T	C, T	C, T	C, T	C, T
RS1821380	C, G	C, G	C, G	C, G	C, G	C, G	C, G	C, G	C, G	C, G	C, G	C, G
RS2073383	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---
RS214955	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G
RS2272998	C, G	C, G	C, G	C, G	C, G	C, G	C, G	C, G	C, G	C, G	C, G	C, G
RS2503107	A, C	A, C	A, C	A, C	A, C	A, C	A, C	A, C	A, C	A, C	A, C	A, C
RS2567608	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G
RS279844	A, ---	A, ---	A, ---	A, ---	A, ---	A, ---	A, ---	A, ---	A, ---	A, ---	A, ---	A, ---
RS315791	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---
RS321198	C, T	C, T	C, T	C, T	C, T	C, T	C, T	C, T	C, T	C, T	C, T	C, T
RS338882	T, ---	T, ---	T, ---	T, ---	T, ---	T, ---	T, ---	T, ---	T, ---	T, ---	T, ---	T, ---
RS3780962	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---	C, ---
RS445251	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---
RS447818	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G
RS560681	A, ---	A, ---	A, ---	A, ---	A, ---	A, ---	A, ---	A, ---	A, ---	A, ---	A, ---	A, ---
RS6444724	C, T	C, T	C, T	C, T	C, T	C, T	C, T	C, T	C, T	C, T	C, T	C, T
RS6591147	C, T	C, T	C, T	C, T	C, T	C, T	C, T	C, T	C, T	C, T	C, T	C, T
RS6811238	T, ---	T, ---	T, ---	T, ---	T, ---	T, ---	T, ---	T, ---	T, ---	T, ---	T, ---	T, ---
RS7205345	C, G	C, G	C, G	C, G	C, G	C, G	C, G	C, G	C, G	C, G	C, G	C, G
RS7229946	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---
RS740598	A, ---	A, ---	A, ---	A, ---	A, ---	A, ---	A, ---	A, ---	A, ---	A, ---	A, ---	A, ---
RS7520386	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G
RS7704770	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---	G, ---
RS985492	T, ---	T, ---	T, ---	T, ---	T, ---	T, ---	T, ---	T, ---	T, ---	T, ---	T, ---	T, ---
RS987640	T, ---	T, ---	T, ---	T, ---	T, ---	T, ---	T, ---	T, ---	T, ---	T, ---	T, ---	T, ---
RS9951171	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G	A, G

Small Test Panel: 96 Individuals

- **10 buccal swab DNAs**
- **38 blood samples obtained commercially**
 - 20 samples tested in parallel with 40 individual TaqMan® assays
- **50 blinded DNA samples from UNTHSC (John Planz)**
 - Tested at UNTHSC with AB GenPlex™ Kidd-40 panel
- **Buccal swabs run at a set dilution factor (corresponded to between 300 pg and 1.8 ng per reaction)**
- **Blood samples were run at 500 pg/reaction**
- **Samples were run in duplicate**
- **Each sample gave a consistent profile between duplicates**
- **100% concordance between TaqMan® and Ibis™ assays**
- **Concordance between Ibis™ and GenPlex™ for all loci except one**
 - rs2073383 was 100% concordant between Ibis™ and TaqMan®
 - rs2073383 discordant between GenPlex™ and TaqMan®

Comparison of 86 Samples Run with Ibis™ Kit to Distribution from Kidd Paper


Pair-wise comparisons between 86 individuals –
3655 comparisons



Number of matching loci between 2 individuals


Comparison of 86 Samples Run with Ibis™ Kit to Distribution from Kidd Paper

Pair-wise comparisons between 86 individuals –
3655 comparisons




Comparison of 86 Samples Run with Ibis™ Kit to Distribution from Kidd Paper

Pair-wise comparisons between 86 individuals –
3655 comparisons



Comparison of 86 Samples Run with Ibis™ Kit to Distribution from Kidd Paper

Pair-wise comparisons between 86 individuals –
3655 comparisons



Summary

- **Basic 40-Kidd-SNP assay defined and moving into validation**
- **Limit of sensitivity down to 63 – 125 pg/reaction**
- **Addition of mammal, fungal or bacterial DNA does not appear to interfere with the assay**
- **100% concordance with AB TaqMan[®] assays for all loci**
- **Small panel of samples showed average locus differences comparable to Kidd results**
- **Integrated software developed and in process**

Questions?

Contact Information

Thomas Hall

Ibis Biosciences division of Abbott Molecular

760-603-2375

thall@ibisbio.com

Note:

All images and charts courtesy of Tom Hall, Ph.D. unless otherwise noted.