



Technology Transition Workshop | *Matt Kramer*

# ***Introduction to Post-PCR Cleanup***

# Overview

- **Why post-PCR amplification cleanup?**
- **Enhancing human identity testing**
  - Introduction to QIAGEN® MinElute® post-PCR cleanup technologies
  - MinElute® as a tool to manage throughput and growth
  - MinElute® as a tool to enhance DNA forensic capabilities and to aid in special cases
- **Customer validation data and examples**
  - Custom fit is key to optimal lab efficiency
  - Online support
- **Automation concepts (tomorrow)**



# ***Why Post-PCR Amplification Cleanup?***

- **Increasing number of touch samples (B & E cases, etc.)**
  - Increased samples --> managing throughput
  - Triaging samples (high/low template DNA)
  - Low copy number (LCN) DNA
    - < 100 pg DNA
    - Few cells, low quantity
- **Increasing demand/need to process specialized, challenging sample types (bone, hair fragments)**
- **Degraded DNA**

# ***Why Post-PCR Amplification Cleanup?***

- **Various techniques to handle special DNA applications**
  - **LCN PCR (increased cycle number and nested PCR)**
  - **MiniFiler™**
  - **Whole genome amplification (WGA)**
  - **Post-amp cleanup**
  - **Others**
    - **mtDNA**
    - **Y-STRs**

# ***Why Post-PCR Amplification Cleanup?***

- **Possible problems to troubleshoot**
  - Allele drop-out
  - Increased stutter
  - Heterozygous peak imbalance
- **Implications**
  - (Re)Validation of STR analysis rules
  - Specialized personnel or routine task, validation and competency testing
  - Ensuring enough reagent blank for the procedure



## Enhancing Human Identity Testing


# *Introduction to QIAGEN<sup>®</sup> Post-PCR Cleanup Technologies*

Technology  
Transition Workshop



# QIAGEN® Kits for Forensic Use

- **Broad range of protocols for forensic casework/ reference samples**
  - **Purification of genomic DNA from...**
    - Hair
    - Cigarette butts
    - Body fluids on swabs
    - Blood, sperm, and saliva stains on fabric
    - Nail clippings and scrapings
    - Chewing gum
    - Urine, feces
    - Tissue biopsy
    - Laser-microdissected tissue
  - **Extraction cleanup of gDNA (P:C)**
  - **Post-PCR amplification cleanup**



# ***Goals of MinElute<sup>®</sup> for Forensic Casework***


- **Increase sensitivity by increasing alleles detected**
  - Enhanced low copy number (LCN) interpretations?
  - Enhanced mixture interpretations?
- **Preserve sample by eliminating need for reduced volume PCR (RV-PCR)**
  - Elimination of increased cycle numbers
- **Reduce complex and expensive techniques like laser-capture microdissection (LCM)**
- **Automation friendly, automate on the QIAcube<sup>®</sup>**







**Unique MinElute  
membrane assembly**



# *Silica Adsorption Nucleic Acid Purification*




# Silica Chemistry – Lysis and Binding




# Wash Away Contaminants With Chaotropic Salt Washes



# Wash Away Chaotropic Salt With Ethanol




# Dry the Silica (Evaporate any Ethanol)



Water molecule

# *Elute in Water or Low Salt Buffer, pH 7-9*



Technology  
Transition Workshop 

JOURNAL OF COLLOID AND INTERFACE SCIENCE **181**, 635–644 (1996)  
ARTICLE NO. 0421

## Driving Forces for DNA Adsorption to Silica in Perchlorate Solutions

KATHRYN A. MELZAK, CHRIS S. SHERWOOD, ROBIN F. B. TURNER,\* AND CHARLES A. HAYNES†<sup>1</sup>

*Biotechnology Laboratory, \*Department of Electrical Engineering, and †Department of Chemical Engineering, 237 Wesbrook Building, 6174 University Boulevard, The University of British Columbia, Vancouver, British Columbia V6T 1Z3, Canada*

- **Shielded intermolecular electrostatic forces**
- **Dehydration of the DNA and silica surfaces**
- **Intermolecular hydrogen bond formation in the DNA-silica contact layer**





# MinElute<sup>®</sup> Kit Contents


## Kit Contents

MinElute PCR Purification Kits	(50)	(250)
Catalog no.	28004	28006
MinElute Spin Columns	50	250
Buffer PB*	30 ml	150 ml
Buffer PE (concentrate)	2 x 6 ml	55 ml
Buffer EB	15 ml	15 ml
pH Indicator	800 µl	800 µl
Collection Tubes (2 ml)	50	250
Loading Dye	110 µl	550 µl
Handbook	1	1

# The MinElute<sup>®</sup> Procedure



# Comparison of QIAGEN® MinElute® Kits




- Purification of > 70 bp amp products
- Primers up to 40 nt long are removed

# MinElute<sup>®</sup> PCR Purification Kit

## Specifications

Product Specifications	
MinElute PCR Purification Kit	
Maximum binding capacity:	5 µg
Recovery of DNA:	80% (70 bp – 4kb)
Maximum weight of gel slice:	—
Elution volume:	10 µl
Volume of eluate:	9 µl
Capacity of column reservoir:	800 µl
<b>Recovered:</b> dsDNA	70 bp – 4 kb
<b>Removed:</b> <40mers	YES

# Mechanism for Post-PCR Purification Enhancement of Trace DNA



# ***Evaluation of Post-PCR Cleanup to Enable LCN Analysis at 28 Amplification Cycles***

- **Smith, Pamela J. and Ballantyne, Jack. “Simplified Low-Copy-Number DNA Analysis by Post-PCR Purification.” *Journal of Forensic Sciences* 52 (4) (July 2007): 820-829.**
  - **Smith and Ballantyne demonstrated for the first time that MinElute<sup>®</sup> post-PCR purification can substantially increase sensitivity, performance**
  - **Simplified LCN DNA analysis where < 100 pg DNA available**
  - **Commonly used techniques to increase sensitivity:**
    - **Increase PCR cycles: 28 to 30-34**
    - **Nested PCR**

# ***Evaluation of Post-PCR Cleanup to Enable LCN Analysis at 28 Amplification Cycles***

- **Smith, Pamela J. and Ballantyne, Jack. “Simplified Low-Copy-Number DNA Analysis by Post-PCR Purification.” *Journal of Forensic Sciences* 52 (4) (July 2007): 820-829.**
  - **Investigate possible issues:**
    - Allele drop-out, drop-in
    - Higher stutter peaks
    - Contamination
  - **Dilution of template DNA from 5 to 600 pg from blood spots on both cotton cloth and FTA cards, analysis in quadruplicates**

# Smith and Ballantyne Study: Comparison of Different Cleanup Procedures

	Signal x Fold	Artifacts (Off-ladder Alleles)	Process Time
Microcon <sup>®</sup> 50	3 to 6	Yes, one locus	55 minutes (4 wash steps)
Montage <sup>®</sup> PCR	6 to 8	Yes, one locus	65 minutes (4 wash steps)
ExoSAP-IT <sup>®</sup>	None	Strong	30 minutes (mostly incubation)
MinElute <sup>®</sup>	4 to 6	None	6 minutes (2 wash steps)



# Smith and Ballantyne Study: Reliable STR Profiles From 20 pg Template Using MinElute®

Same sample analyzed with and without post-amplification purification

STR profile from 20 pg – no post-PCR cleanup




Image courtesy of Smith and Ballantyne

# Smith and Ballantyne Study: Reliable STR Profiles From 20 pg Template Using MinElute®

Same sample analyzed with and without post-amplification purification

Same sample after post-PCR purification and injection of entire concentrated product




Image courtesy of Smith and Ballantyne

# ***Smith and Ballantype Study: Evaluation of Post-PCR Purification LCN Analysis***

- **Low copy number DNA: 5 to 600 pg evaluated**
- **Approach:**
  - **Evaluated effect of MinElute® Post-PCR Purification Kit**
- **Results:**
  - **With MinElute®-based purification full profiles were consistently obtained from 20 pg with 28 PCR cycles**
  - **Significant data obtained down to 5 pg**
  - **Increase in signal intensity 4 to 6-fold (4-fold on average) at constant DNA amount used for injection**
  - **Up to 19-fold signal increase with increased DNA injection amount**
- **Conclusion:**
  - **LCN DNA analysis can be simplified by using post-PCR purification**

## ***Supporting Journal Article***

- **Gill P., Whitaker J., Flaxman C., Brown N., and Buckleton J. “An Investigation of the Rigor of Interpretation Rules for STRs Derived from Less than 100 pg of DNA.” *Forensic Science International* 112 (1) (July 24, 2000): 17-40.**
- **Contact Information for Peter Gill:**
  - **Forensic Science Service, Priory House, Gooch Street North, B56QQ, Birmingham, UK**
  - **[dnapgill@compuserve.com](mailto:dnapgill@compuserve.com)**



**Post-amplification Cleanup  
using the MinElute® Post-PCR Purification Kit:  
*A Presentation of Two More Papers and  
Customer Examples***

# ***Summary of Findings***

- **Post-PCR cleanup can greatly enhance the sensitivity of STR profiling for low copy number casework applications**
- **Improving sensitivity in trace DNA casework applications by post-PCR purification is increasingly receiving attention**
- **Success rates using standard forensic methods to obtain intelligible profiles from trace DNA samples usually range around only 30-50%**
- **Efforts for method enhancements usually focus on the pre-analytical or pre-PCR process of the analysis**
- **However, there can also be tremendous value in post-PCR method optimization**
- **Two publications in leading peer-reviewed journals, in addition to the Smith and Ballantyne paper previously discussed, demonstrate the value of forensic post-PCR purification**

## ***Study 1: Forster, et al.***

- **Forster, L., et al. “Direct Comparison of Post-28-cycle PCR Purification and Modified Capillary Electrophoresis Methods with the 34-cycle “Low Copy Number” (LCN) Method for Analysis of Trace Forensic DNA Samples.” *Forensic Science International - Genetics* 2 (4) (September 2008): 318 – 328.**
- **Forster, et al. challenged increased cycle PCR as a strategy for trace forensic DNA, proposing to combine standard PCR with MinElute® post-PCR purification instead**
- **Direct comparison of post-28-cycle PCR purification and modified capillary electrophoresis methods with the 34-cycle “low copy number” (LCN) method for analysis of trace forensic DNA samples**


# ***Study 1: Forster, et al.***

- **Results:**

- **The combination of 28-cycle PCR with post-PCR purification using the MinElute<sup>®</sup> PCR Purification Kit:**
  - **Provided equivalent results in STR analysis**
  - **Reduced sample consumption**
  - **Reduced stutter peak ratios**
  - **Reduced the frequency of non-donor allele peaks**
- **“This study demonstrates that gains in profile quality at least equivalent to those seen in 34-cycle PCR can be made by concentrating and increasing the loading of product from a 28-cycle PCR.”**



# Study 1: Forster, et al.: 50 $\mu$ l SGM+



**Table 1**  
Shows a summary of the differences between conditions A, B, C and D

Condition	Number of PCR cycles	Volume of sample loaded onto 3130 ( $\mu$ l)	Post-PCR clean-up	Injection time (s)	Injection voltage (kV)
A	28	1	None	10	3
B	28	2	MinElute	10	3
C	28	2	MinElute	30	4
D	34	1	None	10	3

Images courtesy of Forster, et al.

# ***Study 1: Forster, et al.***

- **Summary:**
  - **Trace DNA analysis utilizing standard 28-cycle PCR combined with MinElute<sup>®</sup> post-PCR purification achieves equal or better results than increased cycle number procedures with 34-cycle PCR**
  - **Common problems of increased cycle PCR, such as elevated stutter peak ratios and increased non-donor allele peaks frequencies, can be avoided**

## ***Study 2: Mayntz-Press, et al.***

- **Mayntz-Press, K.A., et al. “Y-STR Profiling in Extended Interval (+/- 3 days) Postcoital Cervicovaginal Samples.” *Journal of Forensic Sciences* 53 (2) (March 2008): 342–348.**
- **Mayntz-Press, et al. applied MinElute<sup>®</sup>-based post-PCR purification after differential extraction in sexual assault analysis**

## ***Study 2: Mayntz-Press, et al.: Results***

- **MinElute<sup>®</sup> purification not only increased peak heights with partial profiles, but also the appearance of loci or allele call rates**
- **Profiles with results at 10 to 17 loci were obtained from 6-day postcoital samples that, prior to purification, yielded only 6 to 11 locus profiles**
- **The incorporation of a simple post-PCR purification process significantly increased the ability to obtain Y-STR profiles, particularly from 5- to 6-day postcoital samples**
- **Post-PCR purification increases the sensitivity of allele detection resulting in:**
  - **Increased peak heights with partial profiles**
  - **The appearance of loci that otherwise were not apparent in the corresponding non-MinElute<sup>®</sup> purified samples**

## ***Study 2: Mayntz-Press, et al.: Study Summary***

- **The post-PCR purification procedure significantly increased the ability to obtain Y-STR profiles**
- **Post-PCR purification increased the sensitivity of allele detection in sexual assault sample analysis**

# ***Study 2: Mayntz-Press, et al.: Continuation at AZ Dept. of Public Safety***

**STRATEGIES FOR INCREASING ALLELE CALLS IN FORENSIC CASEWORK  
USING PCR ENHANCEMENTS AND COMMERCIAL POST AMPLIFICATION  
CLEAN UP SYSTEMS**



Scott C. Milne, B.S., Kathleen A. Mayntz-Press, M.S.F.S.

*Arizona Department of Public Safety, DNA/Serology Casework Unit P.O. Box 6638 MD 1150, Phoenix, AZ 85005*

- **Milne, S.C. and Mayntz-Press, K.A. “Strategies for Increasing Allele Calls in Forensic Casework Using PCR Enhancements and Commercial Post-amplification Cleanup Systems.” (2008) Poster at the 19th International Symposium on Human Identification**

Image courtesy of Milne and Mayntz-Press

Technology  
Transition Workshop 

## ***Study 2: Mayntz-Press, et al.: Continuation at AZ Dept. of Public Safety Results Summary***

- **MinElute<sup>®</sup> PCR Purification Kit results:**
  - Increased the number of loci “called”
  - Showed better sensitivity than other methods compared
  - Useful tool for difficult casework samples
  - Cleaner baseline noted


# ***Study 3: North Louisiana Crime Lab***

- **Study Parameters:**

- **Analysis of pulverized bone powder amplified with the PowerPlex® 16 System kit**
- **Capillary electrophoresis using the Applied Biosystems 3130 Genetic Analyzer**
- **Analysis at 150 RFU using GeneMapper® ID (v.3.2) software**
- **The cleanup procedure was performed as described in the QIAGEN® MinElute® Handbook (03/2008) with a 10 µL elution volume**



# Study 3: North Louisiana Crime Lab



All data generated at the NLCL-Shreveport Headquarters Laboratory

Image courtesy of Jessica Esparza, Ph.D.



Technology  
Transition Workshop

# Study 3: NLCL MinElute<sup>®</sup> Results

## 7.8 pg Template Using PowerPlex<sup>®</sup> 16

Without MinElute<sup>®</sup>

With MinElute<sup>®</sup>




Images courtesy of Jessica Esparza, Ph.D.

Technology  
Transition Workshop




# Study 3: NLCL MinElute® Results

Without MinElute®



With MinElute®



Images courtesy of Jessica Esparza, Ph.D.

## ***Study 3: NLCL MinElute®***

- **Summary of Study Results:**
  - **Total RFU value increased approximately 24-fold after samples subjected to MinElute® cleanup**
  - **Average RFU for alleles detected increased from 241 RFU to 1610 RFU, or approximately 7-fold**
  - **100% of alleles were detected following cleanup compared to only 26% of the alleles detected without MinElute® cleanup**

# Study 4: NFSTC MinElute<sup>®</sup>: Coming Next




Image courtesy of NFSTC

# ***Conclusions: Enhancing Forensic Performance with QIAGEN MinElute Post-PCR Cleanup***

- **Provides ability to enhance DNA forensic performance, data analysis, conclusive results**
- **Provides ability to handle special applications**
  - **Special sample types: bone, hair fragment, etc**
  - **Low template DNA samples (touch, B & E, etc.)**
  - **Degraded DNA**
- **Ability to gain increased RFU values without increased PCR cycles or nested (LV) PCR reactions**
- **Fast, convenient way to improve results, automatable**
- **Comparable or less complex than other techniques (e.g. LCM)**
- **Has to be validated like any other new method**

# Online Support

- **QIAGEN Webinar: Better STR Profiles from Trace Evidence – Validation of MinElute® for Post-amplification Cleanup of Forensic Samples**
  - <http://www1.qiagen.com/events/WebSeminars/MinElute/default.aspx>
  - Jay Caponera of the New York State Police Forensic Investigation Center
    - Validation of the MinElute® procedure
  - Bruce Heidebrecht of Maryland State Police Forensic Sciences Division present on their validation experiences:
    - Validation of the MinElute® procedure on the QIAcube® instrument
- **MinElute® and QIAcube® demo videos can be found at:**
  - <http://www.youtube.com/watch?v=M8TmmL1pcH0>
  - <http://www.youtube.com/watch?v=Li4ISMqSNO8&NR=1>

# *Questions?*



# ***Contact Information***

**Matt Kramer**  
**National Sales Manager**  
**Applied Testing North America**  
**703-623-0564**  
**matt.kramer@qiagen.com**  
**[www.qiagen.com/forensics](http://www.qiagen.com/forensics)**

**Note:** All images are courtesy of QIAGEN® unless otherwise noted.