

Physical Separation of Forensic Mixtures Using

Laser Microdissection Techniques

Technology Transition Workshop

Slide Preparation

1. Pipette 20 μ l of sample cells suspended in **Carnoy's Fixative** onto a membrane slide, 1 mm in thickness, and spread gently into a 1 cm x 1 cm area using a pipette tip.
2. Place the slide above the steam of a water bath for 5 seconds.
3. Place slides onto a slide warmer set to 56°C until completely dry.
4. Bake slides in a biological oven set to 60°C for approximately 3 hours.
5. Store slides in dessicator until ready for pre-treatment.

Slide Pre-Treatment

1. Cover marked sample region with reusable perfusion chamber.
2. Pipette 500 μ l **Pepsin Buffer**, pre-warmed to 37°C, into the chamber. Incubate on the ThermoBrite for 2 minutes at 37°C.
3. Remove the perfusion chamber and place the slide into a Coplin jar containing **PBS/MgCl₂** solution. Incubate at room temperature for 5 minutes. Shake gently.
4. Fix slides in **Formalin Buffer** for 5 minutes at room temperature. Shake gently.
5. Wash the slides by placing them into a Coplin jar containing room temperature **PBS** for 5 minutes.
6. The slides are then dehydrated through a series of **Ethanol Washes** (70%, 85%, and 100%) for one minute each at room temperature. Following dehydration, the slides should be air dried.
7. Store slides in a desiccator until FISH processing is initiated.

Physical Separation of Forensic Mixtures Using

Laser Microdissection Techniques

Technology Transition Workshop

CEP X and Y Probe Hybridization

1. Turn off lights before handling probe.
2. Combine 7 μl of **CEP Hybridization Buffer**[®], 1 μl of **CEP X**[®] **probe**, 1 μl of **CEP Y**[®] **probe**, and 1 μl of **ddH₂O** in a 0.5 μl microcentrifuge tube.
3. Dehydrate the slides through a series of **Ethanol Washes** (70%, 85%, and 100%) for one minute each. These ethanol washes are performed at room temperature. Following dehydration, the slides should be air dried.
4. Apply 10 μl of the **probe mix** directly to the marked sample region and immediately conceal with a cover slip.
5. Seal the cover slip to the slide with rubber cement.
6. Place slide on the ThermoBrite. Run the pre-programmed denaturation/hybridization program for 45 minutes.

Post-Hybridization

1. Following the 45 minutes incubation period, remove slide from the ThermoBrite
2. Using forceps, remove the rubber cement and cover slip.
3. Cover marked sample region with a perfusion chamber and apply 500 μl **SSC/Tween20 Buffer**. Incubate for 2 minutes on the ThermoBrite at 73°C.
4. Remove perfusion chamber. Immerse slide in a Coplin jar containing **SSC Buffer** for one minute at room temperature. Shake gently.
5. Remove slide from Coplin jar. Air dry slide in darkness.
6. Dip slide in a ddH₂O bath 2-3 times.
7. Place 10 μl of diluted **DAPI VectorShield Counterstain** directly onto the marked sample region. Cover with a cover slip. Incubate the specimen in the dark for 15 minutes at room temperature.
8. View the samples using a fluorescence equipped LM microscope with appropriate filters. All three fluorophores may be visualized using a fluorescence capable microscope with DAPI (blue-nucleus), FITC (green-female), and TRITC (red-male) filters.