


A Technique for Microscopical Soil Examinations

Skip Palenik

Microtrace
Elgin, IL USA

www.microtracescientific.com

Microtrace

A close-up photograph of a brown leather boot. A white rectangular tag is affixed to the upper part of the boot. The tag contains the following text: "50001-85 D.L. Mapes" on the first line and "Boots from Suspect" on the second line. The boot shows some wear and a dark stain on the leather. The stitching of the sole is visible at the bottom of the frame.

50001-85 D.L. Mapes
Boots from Suspect


Microtrace

Edmond Locard




Microtrace

Preliminary Separation




Color



Microtrace


Fraction 1





Microtrace

Heavy Mineral Separation



Mineral Separation in Heavy Liquids




Microtrace

Freezing Heavy Minerals in Tip of Tube with LN2



Microtrace

Isolation of Mineral Separates



Heavies Frozen in Bottom of Microcentrifuge Tube



Microtrace

Light Mineral Fraction Washed from Tube




Microtrace

Light and Heavy Fractions Ready for Mounting




Microtrace

PLM Study of Density Fractions




Microtrace

Light Mineral Fraction in 1.540 Cargille Liquid




Microtrace

Plagioclase Feldspars Crossed Polars




Microtrace

Rock Fragments




Microtrace

Heavy Mineral Suite Plane Polarized Light



Microtrace


Mineral Suite Varies by Provenance of Source Rocks



Microtrace



Refractive Index and Dispersion Colors

- Top. Apatite grain in 1.660 refractive index oil. Plane polarized light.
- Bottom. Crossed polars.




Birefringence

- Top. Kyanite in 1.660 index of refraction oil.
- Bottom. Crossed polars





Pleochroism

- Top. Glaucophane
1.660 refractive
index oil. N-S
polarizer.
- Bottom. E-W
polarizer.




Pleochroism

- Top. Euhedral tourmaline in 1.660 index of refraction oil. N-S polars
- Bottom. E-W polars.




Pleochroism

- Top. Subhedral tourmaline in 1.660 refractive index oil.
- Bottom. E-W polars.



Mineral Varieties



- Hypersthene from Mount St. Helens eruption collected from Yakima, WA days after the event. 1.660 index of refraction oil. Top. N-S polars. Bottom. E-W polars.



Microtrace


Mineral Varieties

- Hypersthene from Africa. 1.660 refractive index oil. Top. N-S polars. Bottom. E-W polars.



Mineral Varieties

- Hypersthene from Martinique. 1.660 refractive index oil. Top. E-W polars. Bottom. N-S polars.



Coarse Mineral Fractions

- Examination for identification and surface texture.
- Stereomicroscopy
- Scanning Electron Microscopy
- Cathodoluminescence




Light Microscopy

- Monahan Sand Dunes in West Texas. Mounted in 1.660 refractive index oil for contrast.



Stains and Reagents

- Sahara sand stained with methylene blue to show distribution of amorphous silica (silicic acid) on grain surfaces.



SEM/EDS




Microtrace

SEM of Grain Surfaces

- Rounded quartz grain from Monahan Dunes in Texas showing surface coating.



Indicators on Quartz Grain Surfaces




Microtrace

Diatoms on Marine Quartz Grain Surface




Microtrace

Detail of Diatom on Quartz Grain Surface



Microtrace

“Silica Flowers” Deposited on Quartz Grain Surface



Microtrace

Etching and Dissolution of Silica on Quartz Surface




Microtrace

Deep Etching on Quartz Grain Surface



Microtrace

Fresh Quartz Grain from Glacier in Canada




Microtrace


Cathodoluminescence




Microtrace




Calcite



Wollastonite




Zircon




Willemite

Microtrace

Fraction 2



Acetolysis in 1.5 mL Glass Microcentrifuge Tube



Microtrace



Microtrace


Fraction 2

Identification of Isolates

Optical Microscopy

- Pollen and spores
- Resistant plant tissue
- Tire rubber
- Combustion soot
- Diatoms and plant opal

Pollen Fraction after Acetolysis and Staining



Microtrace

Study of Internal Structure at High Magnification





Examination of Exine Sculpturing



Microtrace

Pollen Fraction without Staining

- Top. Charcoal particle.
- Bottom. Difficultly identified plant tissue.



Pollen Fraction without Staining



Microtrace

Opal Phytoliths

- Distinctive morphology.




Plant Opal



Microtrace

Fraction 3 (clays)



X-Ray Diffraction



Microtrace

FTIR Spectroscopy



Microtrace

Murder Scene Outside Disneyland in California




Microtrace

Concrete Block Substituted for Computers



Microtrace

Heavies Isolated from Concrete Block



Microtrace



Microtrace